

Sliding weight baby scale M102300 CE class III approved

ADE



New improved design: Easy to mount baby tray – without any screws or clips. The soft curved new tray will keep the baby safely in place and allows a quick length measurement, when weighing the infant. The simple and quick to use ADE sliding weighing baby scale is handy for day to day weighing of babies in general situations. Robust and economical with an all steel construction and a precision weighing mechanism.

Technical specification:

Capacity	16 kg
Graduation	10 g
Functions	Manual zero setting and tare
Dimensions tray	575 x 275 x 70 mm
Weight approx.	5,8 kg
Dimensions overall	575 x 300 x 200 mm

