

# Electronic baby scale

M111600-01

**ADE**



MZ111600-002 optional cotton sling

Practical lightweight and compact portable weighing scale for midwives. The extra large foldable handle gives a secure grip and allows weighing of new born babies quick and safe. An optional cotton sling\* will hold the baby comfortable in place. In addition a small carrying bag is optionally available.

## Technical specification:

Capacity:	20 kg
Graduation:	20g
Functions:	Tare, auto-hold, automatic zero setting, units switch kg/lbs, low battery indication
Dimensions:	55 x 22 x 100 mm
Weight approx.:	0,2 kg
Power supply:	2 x CR2032 batteries



MZ111600-003 optional carrying bag with zipper

\*Sling design may differ from picture