

# Electronic baby scale

M118600

**ADE**



The gently contoured baby tray gives generous as well as secure space for large and small infants. The fine graduation of 5 g all through the weighing range of 20 kg guarantees together with the auto-hold function precise and fast weighing results. The large backlight LCD display makes the reading easy also in dimmed surroundings. Simply to operate by four clearly marked soft touch buttons.

## Technical specification:

Capacity	20 kg
Graduation	5 g
Functions	Tare, auto-hold, kg & lbs reading
Dimensions	600 x 280 x 135 mm
Weight approx.	3,2 kg
Power supply	Mains and battery operated

