

Electronic baby scale

M119600

ADE



The ADE Bebea baby weighing scale is the perfect scale for home use to monitor weight changes of new born babies and for light use at clinics. The long term memory function together with the BMIC (Breast-Milk-Intake-Calculator) will help to keep track of the weighing results for a long period. The auto-hold and tare functions make this new baby scale complete. With a weight of only 1,6 kg the scale can be carried easily. An optional carrying bag is available. Operated by standard AAA batteries.

Technical specification:

Capacity:	15 kg
Graduation:	10 g
Dimensions:	600 x 350 x 75 mm
Weight approx.:	1,6 kg
Optional extra:	Carrying bag MZ10030
Functions:	Tare, autohold, automatic zero setting, units switch kg/lbs, low battery indication
Special function:	Long term memory function with BMIC – Brest milk intake calculator

