

Mechanical physician scale  
M318300 CE class III approved

**ADE**



- Solid and robust all steel construction with heavy duty cast iron lever system
- Reading from either side at eye level
- Large base with non slip cover
- Telescopic height rod up to 210 cm by 1 mm graduation
- 2 castors for mobile use

**Technical specification:**

Capacity:	200 kg
Graduation:	100 g
Dimensions:	280 x 560 x 1470 mm
Weight;	15 kg
Platform size:	275 x 375 mm
Measuring range:	60 – 210 cm by 1 mm