

PRODUCT INFORMATION



International Range
Elektronik Floor Scale

Electronic wide body platform and bariatric scale

Low level platform scale with a durable handrail for additional stability and support for the patient. Two wheels allow an easy transport of the unit. Height measurement from 600 to 21000 mm is an integrated part of this scale

- Automatic zero setting
- Tare and pre- tare function
- Low battery indication
- Auto switch off
- Hold function
- BMI
- RS232 interface
- Mains adapter and rechargeable battery operated



M319660-02

Technical Specifications

Capacity	300 kg
Graduation	100 g
Dimension (mm)	520 x 520 x 55
Dimension platform	520 x 520 mm