


Column and floor scales



The perfect scale for daily use at Hospitals, Clinics and GPs offices. Strong and sturdy construction with a total capacity of 250 kg and a fine graduation of 50 g. The turning display allows the reading from both sides of the display. Mains (optionally) and battery operated for maximum flexibility. A BMI function is integrated as standard function.

Technical specification:

Capacity:	250 kg	Function: BMI calculation, mother and child weighing, turning display,
Graduation:	50g<50kg>100 g	hold & auto-hold, auto-switch off
Dimensions:	400 x 350 x 44 mm	Power supply: 4 x AA batteries life time approx. 4000 cycles mains optionally
Weight:	3,8 kg	Levelling feet
Main all off switch		Tap on start function