



The fully automatic mother and child weighing function will allow determining the weight of new born and small children fast and effortless in the arms of a nurse, parent or field worker just by starting the scale via the tap on function. A lightweight splash proof carry bag for transport is supplied as standard. A fully sealed on and off switch is placed beneath the scale to avoid unwanted battery use in transport. The light-coloured slip resistant platform protects the surface from unwanted heating up in bright sun light. The default setting will start the scale in the mother and child weighing mode. A manual switch over to standard weighing is possible by pressing just one button.



Function: Mother and child weighing, automatic zero setting, auto-hold function, auto switch off, levelling feet, slip-proof surface

Power supply: 4 x AA batteries life time approx. 4000 cycles.

Tap on start function.

Battery operation: 180sec. / Solar or mains power 7 Minutes

Optional: Solar power unit with rechargeable batteries & mains adapter.



Technical specification:

- Capacity: 250 kg
- Graduation: 100 g
- Accuracy: +/- 100g
- Dimensions: 346 x 398 x 40-60 mm – levelling feet
- Weight: 3,7 kg
- Main all off switch
- Auto-switch off