



The perfect scale for daily use at Hospitals, Clinics and GPs offices. Strong and sturdy construction with a total capacity of 250 kg and a fine graduation of 100 g. The fully automatic mother and child weighing function will allow determining the weight of new born and small children fast and effortless in the arms of a nurse, parent or field worker just by starting the scale via the tap on function. The display can be turned at a touch of a button towards a helper or the user standing on the scale. Mains (optionally) and battery operated.

**Technical specification:**

Capacity:	250 kg	Function: Mother and child weighing, turning display, automatic zero setting,
Graduation:	100 g	auto switch off
Dimensions:	400 x 350 x 44 mm	Power supply: 4 x AA batteries life time approx. 4000 cycles - mains optionally.
Weight:	3,8 kg	Levelling feet
Main all off switch		Tap on start function