

PRODUCT INFORMATION



International Range Electronic Chair Scales



M400660

Chair Scale

- Large LCD, 25 mm
- Weight unit kg/lbs selectable
- Automatic zero setting
- Tare and hold function
- BMI (Body Mass Index) function
- Power supply: mains and rechargeable batteries
- Arm- and footrests swivelling
- Double brake system

Option: 4 steering wheels M400660-01



M400660-01 – with 4 steering wheels

Technical Specifications

Capacity	250 kg
Graduation	100 g
Dimensions	570 x 935 x 1,055 mm
Weight	approx. 28 kg