


Special purpose scale
PWN5 Compact and nappy scale

ADE



Diaper and nappy scale with a hygienic stainless platform and slightly raised sides. Water resistant according to IP67. The tare function allows subtracting the weight of small containers. Operated by an integrated set of rechargeable batteries and mains.

Technical specification:

Capacity: 5 kg
Graduation: 1 g
Dimensions: 244 x 187 x 54 mm
Weight: 1,8 kg

Function: Tare, automatic zero setting, low battery indication, selectable auto switch off or permanent operation.